

Size = 158X150 mm, Black Print

INSTRUCTIONS FOR USE

VIGGOWAY™
3 Way Stop Cock With Extension Tube

MATERIALS USED:

-Polycarbonate (PC), PP, PVC, HDPE

INDICATIONS:

- To deliver simultaneous infusion of two different IV fluids.
- To administer bolus IV injections through existing access site for IV infusion.
- To take blood samples through existing IV access site, without a second venipuncture.
- For needle-free access in high-risk patients.
- To separate the injection site from infusion site and prevent mechanical irritation of IV access (in patients at high risk for thrombophlebitis).
- To administer IV contrast agents/fluids in patients undergoing CT/MRI/PET scans, or other procedures that require patient to enter a gantry.

CONTRAINDICATIONS:

- Product should not be used in patients with known hypersensitivity to any of the materials used.
- Administration of highly viscous fluids.
- NOT to be used for blood transfusion.

INSTRUCTIONS FOR USE:

ATTACHMENT OF THREE-WAY STOP COCK with EXTENSION TUBE to IV ACCESS

- Ensure that venous access (e.g. IV cannula) is patent by flushing with saline.
- Inspect packing of 3-way Stop Cock with Extension Tube visually to ascertain that package is intact.
- Remove the device from sterile packing taking care to maintain sterility.
- Connect one or two IV Infusion lines to the stop cock ports as required.

Priming:

- Open clamps of IV infusion lines and displace air entirely through ENTIRE LENGTH of EXTENSION TUBE by rotating the tap to open position (all 3 flow directions should be open). If required the priming can be done with normal saline separately to prevent wastage of infusion fluid.

- Rotate the tap to closed position (output from stop cock is closed).
- Remove protector from male luer at the end of extension tube away from stop cock.
- Connect extension tube to venous access device (e.g. IV cannula) and secure the connection by rotating male luer lock.
- Rotate the tap to open position and adjust the flow of infusion fluid.

REMOVAL OF 3-WAY STOP-COCK with EXTENSION TUBE

- Rotate the tap to closed position such that connection with IV access is terminated.
- Clamp input IV infusion lines.
- Keep injection stopper ready for closing IV access (e.g. IV cannula) after removal.
- Unscrew the male luer lock and remove extension tube with stop cock; close IV access with injection stopper or other device.
- Discard the used stop cock and extension tube according to your facility's standard practice.

TERM OF USE :

- Recommended Maximum Duration of Use: Depends on type of medication and other factors; refer to your facility's guidelines or CDC Guidelines:

WARNINGS

The use of this product is restricted to a qualified doctor or a paramedic.
Read instructions before use.
The product should be used according to the instructions for use.

VIGGO MEDICAL DEVICES LIMITED DISCLAIMS ANY RESPONSIBILITY FOR POSSIBLE CONSEQUENCES RESULTING FROM IMPROPER USE.

The product should not be reprocessed.
Visually inspect and carefully check the product and packaging before use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
The product is guaranteed sterile & non-pyrogenic if the package has not been opened or damaged.
Do not clean or re-sterilise. FOR SINGLE USE ONLY. Discard after use.
Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or infections.

Disposal/Discard: Dispose of/Discard the used Device in accordance with your Country's Healthcare and Safety Regulations.

REF

Product Ref. No.

LOT

Batch Number



Date Of Manufacturing



Use By



Do Not Reuse

STERILE EO

Sterilised by Ethylene Oxide Gas



Caution, Consult
Accompanying Documents



Consult Instructions
For Use



Non-Pyrogenic



Do Not Use if
Package is Damaged



Latex free



Do not re-sterilize



Keep away
from sunlight



Rev. Date: 27.02.2014 HR - Drugs - 625-L-B(H)
AW/IFU_AB, Rev. 01